iQ-Check *Listeria* *monocytogenes* II Kit

User Guide

Test for the real-time PCR detection of *L. monocytogenes* in food and environmental samples

Catalog #3578124

BIO-RAD

Table of Contents

Section 1	Introduction.....	1
Section 2	iQ-Check <i>L. monocytogenes</i> II Technology.....	1
Section 3	Kit Components.....	2
Section 4	Shelf Life and Storage.....	2
Section 5	Materials Required but Not Supplied.....	2
	Equipment.....	2
	Supplies.....	3
Section 6	Safety Precautions and Recommendations for Best Results.....	4
Section 7	Protocol.....	5
	Sample Enrichment.....	5
	Free DNA Removal Treatment.....	7
	DNA Extraction.....	7
	Real-Time PCR.....	8
	Data Analysis.....	8
Section 8	Confirmation of Positive Results.....	10
Section 9	Confirmation of Single Colonies Using the iQ-Check Kit.....	10
Section 10	Test Performance and Validations.....	11
Section 11	References.....	12
Section 12	Revision History.....	13
	Appendix — PCR Mix Calculation Guide.....	14

Section 1 Introduction

Conventional bacteriological methods are often long and tedious. In comparison, iQ-Check *Listeria monocytogenes* II is a simple and rapid qualitative test, allowing the detection of specific DNA sequences unique to *L. monocytogenes* found in environmental samples and food products. Using real-time polymerase chain reaction (PCR), *L. monocytogenes*-specific DNA sequences are amplified and detected simultaneously by means of fluorescent probes. Up to 94 samples can be processed, with a minimized risk of contamination and an easy-to-use procedure. The intended users of this kit are trained laboratory personnel who are performing tests to detect *L. monocytogenes*. The use of this test allows results to be obtained within a few hours after enrichment.

Section 2 iQ-Check *L. monocytogenes* II Technology


The iQ-Check *L. monocytogenes* II Kit is a test based on gene amplification and detection by real-time PCR. The kit's ready-to-use PCR reagents contain oligonucleotides (primers and probes) specific for *L. monocytogenes*, as well as DNA polymerase and nucleotides. Detection and data analysis are optimized for use with a Bio-Rad real-time PCR instrument, such as the CFX96 Touch Deep Well Real-Time PCR Detection System.

PCR is a powerful technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling allow DNA denaturation by heat followed by primers annealing to the target region. The DNA polymerase then uses these primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons.

In real-time PCR, specific probes are used to detect DNA during the amplification step by hybridizing to the amplicons. These probes are linked to a fluorophore that fluoresces only when hybridized to the target sequence. FAM is the fluorophore linked to the probe hybridizing to the *L. monocytogenes*-specific DNA sequence. In the absence of target DNA, no fluorescence will be detected. As the amount of amplicons increases with each round of amplification, fluorescence intensity also increases. The optical module measures this fluorescence at the annealing step during each PCR cycle while the associated software plots the fluorescence intensity versus number of cycles.

A synthetic DNA internal control is included in the reaction mix to validate any possible negative results. This control is amplified with a specific probe at the same time as the *L. monocytogenes* target DNA sequence and is detected by a third fluorophore.

This test allows the qualitative detection of *L. monocytogenes* in select foods and environmental samples previously enriched by culture. It includes the following five main steps:

*Please refer to the user guide of the iQ-Check Free DNA Removal Solution (#10000058391) for the conditions of use.

Section 3

Kit Components

The iQ-Check *Listeria monocytogenes* II Kit contains sufficient reagents for 96 tests (94 samples).

Reagent ID	Reagent	Quantity Provided, ml
A	Lysis reagent	1 bottle, 20
B	Fluorescent probes	1 tube, 0.55
C	Amplification mix	1 tube, 4.4
D	PCR negative control	1 tube, 0.5
E	PCR positive control	1 tube, 0.25
F	Lysis beads	1 bottle, 17.6 g

Section 4

Shelf Life and Storage

Once received, the kit must be stored at 2–8°C. Reagents stored at this temperature can be used until the expiration date indicated on the tubes. The shelf life of the lysis reagent is 6 months once it is mixed with lysis beads.

Section 5

Materials Required but Not Supplied

Equipment

- Lab paddle blender for homogenizing test samples
- Incubator for sample microbiological enrichment
- Specific for extraction in sterile 1.5 ml conical screwcap tubes:
 - Benchtop centrifuge (capable of 10,000–12,000 x g)
 - Dry heat block at 37 ± 2°C and/or 95–100°C
- Cell disruptor, such as a Disruptor Genie (Scientific Industries, Inc.)
- Specific for extraction in deep well plate:
 - Heating thermoshaker* capable of maintaining 37 ± 2°C and/or 95–100°C, with a mixing speed of at least 1,300 rpm
- Vortexer
- Magnetic stir plate
- 20, 200, and 1,000 µl micropipets
- Bio-Rad real-time PCR system,* for example, the CFX96 Touch Deep Well System (catalog #3600037) or CFX Opus Deepwell (catalog # 17007991) Real-Time PCR Systems
- Bio-Rad iQ-Check Prep System for automated DNA extraction and PCR plate setup (catalog #3594911)

Note: We recommend using an uninterrupted power supply (UPS) with the thermal cycler and iQ-Check Prep System.

* Contact Bio-Rad Technical Support for information about recommended instruments.

Supplies

- Enrichment medium: *Listeria* Special Broth (LSB) (catalog #3555703, 225 ml x 6 bottles; 3564703, dehydrated, 500 g; 3555793, 5 L x 2 bags; 3564753, dehydrated, 5 kg) or Half Fraser Broth (catalog #3564604, dehydrated, 500 g; 3555797, 225 ml x 6 bottles; 3555794, 3 L x 4 bags)
- Enrichment medium: *Listeria* Special Broth II (LSB II) (catalog #12017463, 225 ml x 6 bottles; 12017388, dehydrated, 500 g; 12017378, dehydrated, 5 kg)
- iQ-Check Free DNA Removal Solution (catalog #3594970)
- Specific for environmental samples:
 - Environmental sponges
 - Environmental swabs
 - Neutralizing broth for sponges and swabs, such as Dey-Engley (D/E) or HiCap Neutralizing Broth, or Lethen Broth
- Specific for extraction in tubes:
 - 1.5 ml conical screwcap tubes, sterile (for example, catalog #2240110XTU)
- Specific for extraction in a deep well plate:
 - 96-well deep well plate (catalog #3594900)
 - Plastic sealing film (catalog #3590139)
 - X-Pierce Sealing Films (catalog #3593977) or Pre-Pierced Plate Sealing Film (catalog #3600040, North America only)
- Specific for Standard II and Easy II extraction protocols:
 - iQ-Check Lysis Beads (reagent F) (catalog #3578136)
 - 200 µl wide-opening pipet tips
- Specific for iQ-Check Prep System:
 - 60 ml dilution container (catalog #3594904)
 - Filter tips (catalog #3594902 or 12014486: 50 µl; 3594903 or 12014483: 1,000 µl)
- PCR mix tubes (catalog #12016673, 5 ml x 25)
- RAPID'*Listeria* spp. Agar* (catalog #3564744, dehydrated, 500 g; 3564745, supplement 1; 3564746, supplement 2)
- RAPID'*L.mono* Agar* (catalog #3563694, 90 mm x 20 dishes; 3555294, kit for 190 ml agar plus 2 supplements; 3564293, dehydrated, 500 g; 3564294, supplement 1; 3564746, supplement 2)
- PCR plates, tubes, sealing tape, and caps
- Sterile filter tips adaptable to 20, 200, and 1,000 µl micropipets
- Tips for Combitip Pipets or equivalent repeat pipettors; sterile, individually packaged
- 1 and 10 ml pipets
- 2 and 5 ml sterile test tubes
- Powder-free gloves
- Distilled sterile water
- Bleach, 5%

- Cleaning agent such as DNA AWAY or RNase AWAY

* Can be used outside of the AOAC-PTM validation.

Section 6

Safety Precautions and Recommendations for Best Results

- This test must be performed by trained personnel
- Pregnant women, children, the elderly, and immunocompromised individuals should neither perform this method nor handle *L. monocytogenes* due to the high infection and fatality rate associated with these groups
- Samples and enrichment cultures must be handled as potentially infectious material and discarded according to local rules and regulations
- All potentially infectious material should be autoclaved before disposal
- The quality of results depends on strict compliance with Good Laboratory Practices (for example, the EN ISO 7218 standard), especially concerning PCR:
 - Never circulate laboratory equipment (pipets, tubes, etc.) from one workstation to another
 - Always use a positive control and a negative control for each series of amplification reactions
 - Do not use reagents after their expiration date
 - Vortex reagents from the kit before using them to ensure homogeneity
 - Periodically verify the accuracy and precision of pipets, as well as the correct functioning of the instruments
 - Change gloves often, especially if you suspect they are contaminated
 - Clean work spaces periodically with 5% bleach and other decontaminating agents, such as DNA AWAY
 - Use powder-free gloves and avoid fingerprints and writing on tube caps. Both will interfere with data acquisition
- It is strongly advised to follow the general requirements described in the standard EN ISO 22174:2005 (Microbiology of food and animal feeding stuffs — Polymerase chain reaction (PCR) for the detection of food-borne pathogens — General requirements and definitions)
- iQ-Check *Listeria monocytogenes* II Kit:
 - All substances or mixtures in the test kit are classified products, according to the Globally Harmonized System (GHS). Contact with acids may cause release of toxic gases. No special precautions are necessary if used correctly. If the product is inhaled, supply fresh air and consult a doctor in case of complaints. After eye contact with the product, rinse opened eye for several minutes under running water. If the products are swallowed, induce vomiting and call for medical help
- iQ-Check Prep System:
 - Improper use of the iQ-Check Prep System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the iQ-Check Prep System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of instrument must be performed only by Bio-Rad field service engineers
- CFX96 Touch Deep Well or CFX Opus Deepwell Real-Time PCR Detection Systems:

- Improper use of the CFX96 Touch Deep Well System or CFX Opus Deepwell System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the CFX96 Touch Deep Well System and CFX Opus Deepwell System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of instrument must be performed only by Bio-Rad Field Service Engineers
- Enrichment:
 - The user should read, understand, and follow all safety information in the instructions for the iQ-Check *Listeria monocytogenes* II Kit. Retain the safety instructions for future reference. To reduce the risks associated with exposure to chemicals and biohazards, perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples. Avoid contact with the contents of the enrichment media and reagent tubes after amplification. Dispose of enriched samples according to current industry standards
 - *Listeria* is a Biosafety Level 2 organism. Biological samples such as enrichments have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations on disposal of biological waste. Wear appropriate protective equipment, which includes, but is not limited to, protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities utilizing the appropriate safety equipment (for example, physical containment devices). Individuals should be trained in accordance with applicable regulatory and company/institution requirements before working with potentially infectious materials
 - When testing is complete, all materials and media possibly containing pathogens should be decontaminated following current industry standards for the disposal of contaminated waste (that is, autoclave for 20 min at 121°C). Consult the Safety Data Sheet for additional information and local regulations for disposal

Section 7 Protocol

A. Sample Enrichment

It is strongly recommended to read the entire protocol before starting the test.

Enrichment media must be at the appropriate incubation temperature (ambient or 30°C/37°C when required) before use.

The use of an enrichment bag with incorporated filter is highly recommended.

Within the scope of the NF VALIDATION, test portions weighing more than 25 g have not been tested. Incubate, without shaking, for times and at temperatures indicated in the table below. It is also possible to carry out the iQ-Check test from LSB that has been enriched and then stored at 2–8°C for up to 72 hr.

The following table outlines the different protocols that can be used, depending on the application and the scope of the validation:

NF VALIDATION BRD 7/10-04/05 and NordVal 037		
Scope (Matrices)	Sample Preparation	Enrichment/DNA Extraction
All food products and environmental samples ¹	Homogenize <i>n</i> g food in 9 x <i>n</i> ml Half Fraser Broth or LSB	Standard II extraction (Half Fraser) <ul style="list-style-type: none"> Incubate 24 – 26 hr at 30 ± 1°C Tube format Standard II extraction (LSB) <ul style="list-style-type: none"> Incubate 22–24 hr at 30 ± 1°C Tube format Easy II extraction (LSB) <ul style="list-style-type: none"> Incubate 24–26 hr at 30 ± 1°C Tube format
Production environment and surfaces ²	Homogenize <i>n</i> g food in 9 x <i>n</i> ml prewarmed LSB or room temperature LSB II	Easy II extraction (LSB) <ul style="list-style-type: none"> Incubate 18–26 hr at 30 ± 1°C Tube/deep well format Easy II extraction (LSB II) <ul style="list-style-type: none"> Incubate 18–26 hr at 37 ± 1°C Tube/deep well format
Composite food ²	Homogenize <i>n</i> g food in 9 x <i>n</i> ml LSB II	Easy II extraction (LSB II) <ul style="list-style-type: none"> Incubate 18-26 hr at 37 ± 1°C Deep well format
Dairy products ²	Homogenize <i>n</i> g food in 9 x <i>n</i> ml prewarmed LSB II	Easy II extraction (LSB II) <ul style="list-style-type: none"> Incubate 20-28 hr at 37 ± 1°C Deep well format
AOAC PTM 010802		
Scope (Matrices) ^{2,3}	Sample Preparation	Enrichment/DNA Extraction
Deli turkey, hot dogs, smoked salmon, cottage cheese	Homogenize <i>n</i> g food in 9 x <i>n</i> ml LSB or LSB II (25 g in 225 ml; 125 g in 1,125 ml)	Easy II extraction (LSB) <ul style="list-style-type: none"> Incubate 24–26 hr at 30 ± 1°C Tube/deep well format Easy II extraction (LSB II) <ul style="list-style-type: none"> Incubate 18-24 hr at 37 ± 1°C Tube/deep well format
Stainless steel surfaces	Environmental samples <ul style="list-style-type: none"> Moisten swabs and sponges with a neutralizing broth that does not contain aryl sulfonate complex For surfaces being analyzed with swabs, sample a 1 x 1" (2.54 x 2.54 cm) area For surfaces being analyzed with sponges, sample a 4 x 4" (10.16 x 10.16 cm) area Add to enough LSB II to cover the swab (10 ml) or sponge (60-225 ml) 	Easy II extraction <ul style="list-style-type: none"> Incubate 16-24 hr at 37 ± 1°C Tube/deep well format
Refer to MFLP-54 for Health Canada approved matrices and protocols		

¹ Validation includes “Lmono Fast” Application Protocol File for a reduced PCR run time. Contact your Bio-Rad representative for more information.

² Validation includes the use of the iQ-Check Free DNA Removal Solution and of the “Lmono Fast” Application Protocol File for a reduced PCR run time. Contact your Bio-Rad representative for more information.

³ Validation includes direct streak to RAPID *L. mono*, RAPID *Listeria* spp., or Agar *Listeria* plates (excluding natural cheese).

B. Free DNA Removal Treatment

The iQ-Check Free DNA Removal Solution provides an ideal way to remove free DNA.

Follow Bio-Rad's recommendations in the user guide (document #10000058391). This protocol is validated by AFNOR Certification in the production environmental samples category.

C. DNA Extraction

General recommendations:

1. Turn on the heat block or thermoshaker to preheat before starting the test. Set it to 95–100°C. Keep the lysis reagent in suspension while pipetting by stirring at medium speed on a magnetic stir plate.
2. In general, avoid shaking the enrichment bag and collecting large fragments of food debris. For food samples with a fatty supernatant, collect the sample just below this layer.
3. Open tubes and wells carefully to avoid any possible cross-contamination.
4. Cool the deep well plate before pipetting directly through pre-pierced sealing film.
5. Use the magnetic bar to keep the lysis reagent in suspension. Pipet while it is stirring at medium speed.
6. Gently shake the lysis reagent by hand first to resuspend the resin. Then pipet while it is stirring at medium speed with the magnetic bar contained in the bottle, in order to keep it in suspension.
7. Reconstitute the final lysis reagent as follows:
 - a. Carefully pour all the contents from reagent F (lysis beads) into reagent A (lysis reagent).
 - b. Use consumables with a wide enough tip to allow pipetting of the homogenized lysis reagent.
 - c. Lysis reagent mixed with lysis beads (reagents A + F) has a shelf life of 6 months when stored at 4°C.

Standard II Protocol

1. Collect 1.5 ml of decanted enriched sample into a tube.
2. Centrifuge the tube at 10,000–12,000 x g for 5 min. Discard the supernatant.
3. Add 250 µl of homogenized lysis reagent (reagents A + F) to the pellet and resuspend the pellet by pipetting the reagent up and down. Close the tube.

Note: Gently shake the lysis reagent by hand first to resuspend the beads.

4. Place the tube in the thermoshaker or cell disruptor for 3 ± 1 min.
5. Incubate in the appropriate heat block at 95–100°C for 10–15 min.
6. Vortex the tube at high speed and centrifuge at 10,000–12,000 x g for 5 min.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

Easy II Protocol

1. Aliquot 100 µl of homogenized lysis reagent (reagents A + F) to the tubes or wells of a deep well plate.

Note: Gently shake the lysis reagent by hand first to resuspend the beads.

2. Add 100 µl of enriched sample.

Note: Shake the suspension to homogenize the culture and then allow any debris to settle before collecting the sample.

3. Mix the solution by pipetting up and down until homogenized.
4. Close the tubes or seal the deep well plate with pre-pierced sealing film.
5. Place the tubes in the cell disruptor for 3 ± 1 min. Incubate the tubes in the heat block at $95\text{--}100^\circ\text{C}$ for 15–20 min.
6. If using the agitator-incubator for the tube or deep well plate format, incubate the tubes or deep well plate in the agitator-incubator at 1,300 to 1,600 rpm at $95\text{--}100^\circ\text{C}$ for 15–20 min.
7. Vortex the tubes at high speed, and then centrifuge at $10,000\text{--}12,000 \times g$ for at least 2 min.

Centrifugation is not needed for the deep well plate.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at -20°C . Always allow it to thaw and homogenize, and then centrifuge at $10,000\text{--}12,000 \times g$ for 5 min before reusing.

D. Real-Time PCR

Instrument and Software Setup

For instrument and software setup, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

PCR Mix Preparation

1. Prepare the PCR mix containing the amplification solution (reagent C) and the fluorescent probes (reagent B). The volume of PCR mix needed depends on the number of samples and controls to be analyzed. At least one positive and one negative control must be included in each PCR run. Use the pipetting table in Appendix — PCR Mix Calculation Guide to find the correct volumes to use for each reagent.

Note: Use the PCR mix (reagents B + C) immediately after preparation. It is stable for 1 hr maximum at $2\text{--}8^\circ\text{C}$.

2. Pipet $45 \mu\text{l}$ of the PCR mix into each well according to your plate setup.
3. Add $5 \mu\text{l}$ of DNA extract, reagent D (negative control), or reagent E (positive control). Do not vortex the sample before pipetting. Hermetically seal the wells of the PCR plate or tube strips. It is important to avoid bubbles at the bottom of the wells by pipetting carefully. As an optional step, centrifuge the sealed PCR plate or tube strips (quick spin) to eliminate any bubbles.
4. Place the PCR plate or tube strips in the thermal cycler. Be sure to place the plate correctly, with the A1 well at the upper left corner. Close the reaction module.

Run PCR

To start the PCR run, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

The “Lmono Fast” Application Protocol File has been NF and AOAC validated.

E. Data Analysis

Data can be analyzed directly at the end of the PCR run or at a later time by opening the stored data file. Follow instructions in the corresponding CFX Manager IDE Software User Manual for opening data files and setting the data analysis parameters.

Interpreting Results

Once the data analysis parameters have been set, results are interpreted by analyzing the quantification cycle (C_q) values of each sample (the cycle at which the amplification curve crosses the threshold).

CFX Manager IDE Software allows complete automated analysis for Bio-Rad real-time PCR detection systems. A verification of the typical characteristics of the amplification curves should be performed prior to releasing results. Contact your Bio-Rad technical support team if additional support is required.

Controls

Verify the positive and negative controls before interpreting sample results.

For the experiment to be valid, the controls must have the results summarized in the following table. Otherwise the PCR reaction must be repeated.

Control	<i>L. monocytogenes</i> Detection (FAM Channel)	Internal Control Detection (HEX Channel)
Negative control	Cq = N/A*	28 ≤ Cq ≤ 40
Positive control	26 ≤ Cq ≤ 36	N/A

* The software indicates a Cq value (the cycle at which the amplification curve crosses the threshold) of N/A (not applicable) when the fluorescence of a sample does not rise significantly above the background noise, and hence does not cross the threshold.

If the results of negative and positive controls differ from those in the Controls table (invalid control), repeat the run and analysis described in D. Real-Time PCR and E. Data Analysis in Section 7 Protocol

Samples

A **positive** iQ-Check *Listeria monocytogenes* PCR test must show a typical amplification curve and a Cq value ≥10 for the FAM fluorophore.

- If the Cq value for both channels is below 10, verify the raw data curve is a regular amplification curve (with a flat baseline followed by a rapid exponential increase of fluorescence and then a flattening out). If the curve seems correct, it may be considered a positive *Listeria monocytogenes* test

If there is no Cq value (Cq = N/A) for FAM, or if the curve is not a typical amplification curve, the internal control for that sample must then be analyzed:

- This sample is considered a **negative** *Listeria monocytogenes* sample if there is no Cq value in the FAM channel and the internal control has a Cq ≥28
- Should the internal control also not have a Cq value (Cq = N/A), this probably indicates an inhibition of the PCR reaction. The sample needs to be diluted (using 10 µl of DNA extract, perform a 1:10 dilution in distilled sterile water, and then test 5 µl of the dilution) and the PCR repeated
- Should the Cq value for the internal control be <28, it is not possible to interpret the result. Verify that the threshold was correctly placed or that the curve as raw data is a regular amplification curve. If the curve does not have a characteristic shape, it will be necessary to repeat the PCR test.

The interpretation of test results is summarized in the following table.

<i>L. monocytogenes</i> Detection (FAM Channel)	Internal Control Detection (HEX Channel)	Interpretation
Cq ≥ 10	N/A	Positive
Cq = N/A	Cq ≥ 28	Negative
Cq = N/A	Cq = N/A	Inhibition*

* When both target and internal control detection give a Cq value = N/A, the DNA extract must be diluted 1:10 and tested again.

An invalid interpretation can be given when validation criteria are not met. Check the raw data and proceed as if the sample was inhibited.

Section 8

Confirmation of Positive Results

In the context of the NF VALIDATION certified method, all positive iQ-Check results must be confirmed in the following ways:

1. Using classic tests, described in the standardized methods CEN or ISO, directly from the Half Fraser Broth, or starting back from the sample when using LSB or LSB II.
2. By isolation of enriched LSB or LSB II (streaking 100 µl) on RAPID'*Listeria* spp. or RAPID'*L.mono* chromogenic agar and incubation for 24 hr at 37 ± 1°C. The presence of characteristic colonies on RAPID'*L.mono* is sufficient to confirm the presence of *L. monocytogenes*. Use an identification test to confirm *L. monocytogenes* on RAPID'*Listeria* spp.
3. By isolation of enriched LSB or LSB II (streaking 10 µl) on AL or PALCAM agar and incubation for 24–48 hr at 37 ± 1°C. Typical colonies should be confirmed by the tests described in the ISO method until species is identified.
4. Using any other method certified by NF VALIDATION based on a principle different from that used in the iQ-Check *L. monocytogenes* II PCR test. The validated protocol of this second method must be followed entirely.

In the case of discrepant results between the iQ-Check *Listeria monocytogenes* II Kit and any of the confirmation options listed above, follow the necessary steps to ensure valid results.

It is possible to store the enriched LSB or LSB II at 2–8°C for 72 hr maximum, following the incubation at 30°C, before carrying out the confirmation.

In the context of AOAC validation, a positive iQ-Check *L. monocytogenes* result should be considered presumptive positive and confirmed according to an appropriate reference method (for example, USDA MLG, FDA BAM, ISO, MFHPB, etc.) is recommended. Alternatively, for select matrices, streak 10 µl enriched LSB or LSB II directly to RAPID'*L. mono*, RAPID'*Listeria* spp. or Agar *Listeria* chromogenic agar and incubate for 24 hr at 37 ± 1°C. The presence of characteristic *Listeria monocytogenes* colonies should be biochemically confirmed for the presence of *Listeria monocytogenes*.

Section 9

Confirmation of Single Colonies Using the iQ-Check Kit

The iQ-Check *Listeria monocytogenes* II Kit may also be used to confirm single isolated *L. monocytogenes* colonies on agar plates. This is officially validated by AFNOR Certification on RAPID'*L.mono* and RAPID'*Listeria* spp. chromogenic media.

1. Pick an isolated colony from a selective or nonselective agar plate with a toothpick, sterile loop, or other adapted consumable (for example, a pipet tip).
2. Resuspend the colony in 100 µl tryptone salt or distilled sterile water in a microcentrifuge tube.
Homogenize using a vortexer.
3. Use 5 µl of the suspension with 45 µl of PCR mix (see D. Real-Time PCR in Section 7 Protocol) and follow the rest of the iQ-Check *L. monocytogenes* II Protocol for the data and results interpretation.

Section 10

Test Performance and Validations

The iQ-Check *Listeria monocytogenes* II Kit is specific for the detection of *L. monocytogenes*

BRD 07/10-04/05

ALTERNATIVE
ANALYSIS
METHODS FOR
AGRIBUSINESS

<http://nf-validation.afnor.org>

NF Validation

The iQ-Check *Listeria monocytogenes* II Kit is certified NF VALIDATION as an alternative method to the reference method ISO 11290-1:2017 (see reference 4 in Section 11), for the detection of *L. monocytogenes* in all products for human consumption and environmental samples. The validation followed the protocol of ISO 16140-2:2016 (see reference 5 in Section 11) and includes the use of the CFX96 Touch Deep Well and CFX Opus Deepwell Real-Time PCR Systems. The use of the iQ-Check Free DNA Removal Solution is validated for composite foods, dairy products and production environmental samples. The associated software is CFX Manager IDE Software (version 2. 2 and later). The “Lmono Fast” Application Protocol File is validated. Certificate number: BRD 07/10–04/05. Valid until: Refer to the certificate available on the AFNOR Certification website.

AOAC-Validation

The iQ-Check *Listeria monocytogenes* II Kit (Easy II Protocol) is validated by the AOAC Research Institute under the Performance Tested Method Program for the detection of *L. monocytogenes* in deli turkey, hot dogs, smoked salmon, cottage cheese, and stainless steel surfaces. A positive result with iQ-Check should be considered presumptive and it is recommended that it be confirmed following the recommendation in Section 8. LSB, LSB II, the “Lmono Fast” APF, the use of the iQ-Check Free DNA Removal Solution, and the use of CFX96 Touch Deep Well and CFX Opus Deepwell Real-Time PCR Systems are validated for all samples. The associated software is CFX Manager IDE Software (version 2.2 and later). Certificate number: 010802.

Health Canada-Validation

The iQ-Check *Listeria monocytogenes* II Kit (Easy II Protocol) has been validated by Health Canada (MFLP-54) for detection of *L. monocytogenes* from RTE meat and poultry, fruit- and vegetable-based products (excluding raw processed vegetables), fish and seafood products (excluding smoked fish), and frozen and fermented dairy products for 25 g and 125 g sample sizes. A positive result with iQ-Check should be considered presumptive and must be confirmed according to MFHPB-30 (see reference 2 in section 11).

NordVal-Validation

The iQ-Check *Listeria monocytogenes* II Kit is validated by NordVal for:

- Enrichment in Half Fraser Broth followed by the Standard II Lysis Protocol
- Enrichment in *Listeria* Special Broth followed by the Standard II Lysis Protocol
- Enrichment in *Listeria* Special Broth followed by the Easy II Lysis Protocol, no longer requiring the first centrifugation step

Certificate number: 037. Refer to the certificate available on the NordVal website for expiration.

Section 11

References

AOAC Official Method 993. 12-1996(1999). *Listeria monocytogenes* in milk and dairy products. Selective enrichment and isolation method.

Health Canada (2011). Health Products and Food Branch. MFHPB-30. Isolation of *Listeria monocytogenes* and other *Listeria* spp. from foods and environmental samples. The Compendium of Analytical Methods, Volume 2, February 2011.

ISO 7218. Microbiology of food and animal feeding stuffs — General requirements and guidance for microbiological examinations.

ISO 11290-1:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 1: Detection method.

ISO 16140-2:2016. Microbiology of the food chain — Method validation — Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method.

ISO 22174. Microbiology of food and animal feeding stuffs — Polymerase chain reaction (PCR) for the detection of food-borne pathogens — General requirements and definitions.

United States Department of Agriculture, Food Safety and Inspection Service (2021). Microbiology Laboratory Guidebook. Chapter 8. 13: Isolation and identification of *Listeria monocytogenes* from red meat, poultry, ready-to-eat siluriformes (fish) and egg products, and environmental samples.

https://www.fsis.usda.gov/sites/default/files/media_file/2021-09/MLG-8.13.pdf, accessed February 14, 2020.

United States Food and Drug Administration (2017). Bacteriological Analytical Manual. Chapter 10: Detection of *Listeria monocytogenes* in foods and environmental samples, and enumeration of *Listeria monocytogenes* in foods. <https://www.fda.gov/food/laboratory-methods-food/bam-chapter-10-detection-listeria-monocytogenes-foods-and-environmental-samples-and-enumeration>, accessed February 14, 2020.

Section 12

Revision History

Release date	Document number	Change
October 2023	10000167776 Ver A	<ul style="list-style-type: none">- Document number change - previous version 10000123535 Ver C iQ-Check <i>L. monocytogenes</i> User Guide- AFNOR Extension: LSB II for composite food, environmental production samples and dairy products

Appendix — PCR Mix Calculation Guide

To find the correct volumes for preparing the PCR mix, add the total number of samples and controls to be analyzed and find the corresponding volumes of reagent B and reagent C in the table.

Total Number of Samples and Controls	Probes Reagent B, μ l	Amplification Mix Reagent C, μ l	Total Number of Samples and Controls	Probes Reagent B, μ l	Amplification Mix Reagent C, μ l	Total Number of Samples and Controls	Probes Reagent B, μ l	Amplification Mix Reagent C, μ l
1	5	40	33	178	1,400	65	351	2,800
2	11	86	34	184	1,500	66	356	2,900
3	16	130	35	189	1,500	67	362	2,900
4	22	173	36	194	1,600	68	367	2,900
5	27	216	37	200	1,600	69	373	3,000
6	32	259	38	205	1,600	70	378	3,000
7	38	302	39	211	1,700	71	383	3,100
8	43	346	40	216	1,700	72	389	3,100
9	49	389	41	221	1,800	73	394	3,200
10	54	432	42	227	1,800	74	400	3,200
11	59	475	43	232	1,900	75	405	3,200
12	65	518	44	238	1,900	76	410	3,300
13	70	562	45	243	1,900	77	416	3,300
14	76	605	46	248	2,000	78	421	3,400
15	81	648	47	254	2,000	79	427	3,400
16	86	691	48	259	2,100	80	432	3,500
17	92	734	49	265	2,100	81	437	3,500
18	97	778	50	270	2,200	82	443	3,500
19	103	821	51	275	2,200	83	448	3,600
20	108	864	52	281	2,200	84	454	3,600
21	113	907	53	286	2,300	85	459	3,700
22	119	950	54	292	2,300	86	464	3,700
23	124	994	55	297	2,400	87	470	3,800
24	130	1,000	56	302	2,400	88	475	3,800
25	135	1,100	57	308	2,500	89	481	3,800
26	140	1,100	58	313	2,500	90	486	3,900
27	146	1,200	59	319	2,500	91	491	3,900
28	151	1,200	60	324	2,600	92	497	4,000
29	157	1,300	61	329	2,600	93	502	4,000
30	162	1,300	62	335	2,700	94	508	4,100
31	167	1,300	63	340	2,700	95	513	4,100
32	173	1,400	64	346	2,800	96	518	4,100

Visit [bio-rad.com/iqcheck](https://www.bio-rad.com/iqcheck) for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

IQ-CHECK is a trademark of Bio-Rad Europe GmbH in certain jurisdictions.

All trademarks used herein are the property of their respective owners.

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United States Department of Agriculture, Food Safety and Inspection Service (2021). Microbiology Laboratory Guidebook. Chapter 8 .13: Isolation and identification of *Listeria monocytogenes* from red meat, poultry, ready-to-eat siluriformes (fish) and egg products, and environmental samples.
https://www.fsis.usda.gov/sites/default/files/media_file/2021-09/MLG-8.13.pdf, accessed February 14, 2020.

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Sezione 11

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AOAC Official Method 993.12-1996(1999). *Listeria monocytogenes* in milk and dairy products. Selective enrichment and isolation method.

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United States Department of Agriculture, Food Safety and Inspection Service (2021). Microbiology Laboratory Guidebook. Chapter 8 .13: Isolation and identification of *Listeria monocytogenes* from red meat, poultry, ready-to-eat siluriformes (fish) and egg products, and environmental samples.

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Seção 11

Referências

AOAC Official Method 993 .12-1996(1999). *Listeria monocytogenes* in milk and dairy products. Selective enrichment and isolation method.

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ISO 11290-1:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Parte 1: Detection method.

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United States Department of Agriculture, Food Safety and Inspection Service (2021). Microbiology Laboratory Guidebook. Chapter 8 .13: Isolation and identification of *Listeria monocytogenes* from red meat, poultry, ready-to-eat siluriformes (fish) and egg products, and environmental samples. https://www.fsis.usda.gov/sites/default/files/media_file/2021-09/MLG-8.13.pdf, accessed February 14, 2020.

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