Health Products and Food Branch

Ottawa

Detection of *Listeria monocytogenes* from selected foods using iQ-Check™ *Listeria monocytogenes* Real-Time PCR Test Kit

Microbiological Methods Committee
Microbiology Evaluation Division
Bureau of Microbial Hazards, Food Directorate,
Health Products and Food Branch, Health Canada
Postal Locator: 2204E
Ottawa, Ontario K1A 0K9

Contact the Microbiological Methods Committee: mmc-cmm@hc-sc.gc.ca

1. Application

This method is applicable to the rapid detection of *Listeria monocytogenes* to determine compliance with the requirements of Sections 4 and 7 of the *Food and Drugs Act* and/or other relevant federal regulations. This method has been validated for use in ready-to-eat meat and poultry, fruit and vegetable based products (except raw processed vegetables), fish and seafood products (except smoked fish), and frozen and fermented dairy products.

2. Description

The iQ-Check™ system is a real-time PCR based assay that involves simultaneous amplification and detection of target organism after a 25 ± 1 hr enrichment in *Listeria* Special Broth. Results are automatically determined by the iQ-Check™ software.

Note: While this method is only approved for certain food products, as listed above, it is assumed that this method could be used with other foods. To ensure the method is fit for purpose for commodities outside the application, it is imperative that other commodities be properly validated following the criteria in the Compendium of Analytical Methods. It is requested that these validation data be forwarded to the [Microbiological Methods Committee](#) (MMC) so the Application Section can be expanded to include these new foods if the data comply with MMC requirements (refer to Development of Methods in [Volume 1 of the Compendium of Analytical Methods](#)).

3. Principle

The iQ-Check™ *Listeria monocytogenes* kit is based on gene amplification and detection by real-time PCR. Ready-to-use PCR reagents contain DNA polymerase, nucleotides, DNA primers and a fluorescent DNA hybridization probe specific to a highly conserved region on the genome of *Listeria monocytogenes*. All reagents necessary for the test are contained in the kit. In real-time PCR, probes are linked to a fluorophore which fluoresces only when hybridized to the target sequence. In the iQ-Check™ *Listeria monocytogenes* kit, carboxyfluorescein (FAM) is the

fluorophore linked to the probe hybridizing to the *Listeria monocytogenes* specific DNA sequence. In the absence of target DNA, no fluorescence will be detected, and the sample is determined to be negative by the iQ-Check™ software. As the amount of amplicon increases with each round of amplification, fluorescence intensity also increases. During each PCR cycle, at the annealing step, the real-time PCR system measures this fluorescence and the associated software plots the fluorescence intensity versus number of cycles. This method allows a simple determination of the presence of *Listeria monocytogenes* in a sample. To monitor for successful DNA amplification in each reaction well, a synthetic DNA "internal control" is included in the reaction mix to monitor for potential inhibition. This control is amplified with a specific probe at the same time as the target DNA sequence, and detected by a second fluorophore.

4. Definition of Terms

See [Appendix A of Volume 1](#).

5. Collection of Samples

See [Appendix B of Volume 1](#).

6. Materials and Special Equipment

Note: The laboratory supervisor must ensure that completion of the analysis described in this method is done in accordance with the International Standards reference "ISO/IEC 17025:2005 (or latest version): General Requirements for the Competence of Testing and Calibration Laboratories".

Note: It is the responsibility of the laboratory to ensure equivalency if any variations of the media formulations listed here are used (either product that is commercially available or made from scratch). Please forward equivalency data to the [Editor of Compendium of Analytical Methods](#) for consideration of modification of this method.

1) iQ-Check™ *Listeria monocytogenes* kit (Bio-Rad catalog # 3578124)

Reagent A - Lysis reagent, 1 bottle (20 mL), contains magnetic stirbar
Reagent B - Fluorescent probes, 1 tube (0.55 mL)
Reagent C - Amplification mix, 1 tube (4.4 mL)
Reagent D - PCR negative control, 1 tube (0.5 mL)
Reagent E - PCR positive control, 1 tube (0.25 mL)
Reagent F - Lysis beads, 1 tube (17.6g)

2) Supplies and reagents

Listeria Special Broth (LSB, Bio-Rad catalog # 3555703 for 225mlx6, # 3564703 for 500g or # 3564753 for 5kg)
Conical screw cap sterile tubes, 1.5 mL (for tube extraction only)
Deepwell microplates (for deepwell extraction only)
Low-profile unskirted PCR plates, 48 or 96 well, white (Bio-Rad catalog# MLL4851XTU, MLL9651XTU, HSP9955XTU)
Optical Flat 8-cap strips, ultraclear (Bio-Rad catalog # TCS0803XTU)
Pre-pierced sealing film (Bio-Rad catalog # 3600040) or equivalent (for deepwell extraction only)
Powder-free gloves

Sterile filter tips - adaptable to 20 μ L, 200 μ L and 1000 μ L micropipettes
Stomacher bag or stomacher bag with incorporated filter (recommended for meat and highly particulate samples)
Micropipettes 20 μ L, 200 μ L and 1000 μ L
PCR mix tube for iQ-Check™ Prep system (Bio-Rad catalog # 3594901)
Pipette tips for iQ-Check™ Prep system (Bio-Rad catalog # 3594902 and 3594903)
Plating and confirmation media and supplies (i.e., from MFHPB-30)
Positive and negative control cultures (e.g., ATCC strains)

3) Equipment

Bench top centrifuge - maximum 10,000-12,000 \times g, for 1.5 mL tubes (for tube extraction only)
CFX Manager Software Industrial Diagnostic Edition (Bio-Rad catalog # 3600042)
Dry heat block, for 1.5 mL tubes - capable of maintaining 95 - 100°C (for tube extraction only)
Incubator - capable of maintaining temperatures of 30 °C for sample incubation
iQ-Check Prep automation system (optional) – for automated DNA extraction and PCR plate setup (Bio-Rad catalog # 3594911)
Magnetic stir plate
Plate sealer – for sealing hard shell PCR plates (Bio-Rad catalog # 1814000)
Real-time PCR system – CFX96 real-time PCR detection system (Bio-Rad catalog # 3600037) or MiniOpticon real-time PCR detection system (Bio-Rad catalog # 3591592)
Stomacher, blender or equivalent
Thermomixer® – capable of maintaining 95-100°C (Bio-Rad catalog # 3600038 [base] and 3600039 [deepwell block] or catalog # 3594999) (for deepwell extraction only)
Vortex mixer

Note: It is the responsibility of each laboratory to ensure that the temperatures of the incubators or water baths are maintained at the recommended temperatures. The following applies to steps of the method which apply to growth only. Where a temperature of $\leq 37^\circ\text{C}$ is recommended in the text of the method the incubator may be at $35 \pm 1.0^\circ\text{C}$. However, where higher temperatures are recommended, it is imperative that the incubators or water baths be maintained within 0.5°C due to potential lethality of the higher temperatures on the microorganism(s) being isolated.

7. Procedure

The test shall be carried out in accordance with the following instructions:

7.1 Handling of Sample Units

7.1.1 In the laboratory prior to analysis, except for shelf-stable foods, keep sample units refrigerated or frozen, depending on the nature of the product. Thaw frozen samples in a refrigerator, or under time and temperature conditions which prevent microbial growth or death.

7.1.2 Analyze sample units as soon as possible after their receipt in the laboratory.

7.2 Preparation of sample

To ensure a representative analytical unit, agitate liquids or free flowing materials until the contents are homogeneous. If the sample unit is a solid, obtain the analytical unit by taking a portion from several locations within the sample unit.

- 7.2.1 Temper LSB to 30°C before use.
- 7.2.2 Prepare a 1 in 10 dilution of the food by aseptically adding 25 g or ml into 225 ml of LSB.
- 7.2.3 Other analytical unit sizes, such as 65 g or ml may be used while maintaining the 1 in 10 dilution. Do not exceed the maximum composited sample size of 125 g or ml.
- 7.2.4 Blend, stomach or vortex as required for thorough mixing.
- 7.2.5 Incubate all samples 30°C for 25 ± 1 hr.

7.3 Preparation for Analysis

- 7.3.1 Turn on heating block (for tube extraction) or Thermomixer[®] (for deepwell extraction) in advance to reach 95 - 100°C.
- 7.3.2 Turn on the thermal cycler and computer.
- 7.3.3 Turn on iQ-Check[™] Prep automation system (for automated sample processing)
- 7.3.4 Carefully pour all the contents from Reagent F (lysis beads) into Reagent A (lysis reagent).
- 7.3.5 Proceed to either step 7.4 for DNA extraction with the tube format or step 7.5 for the DNA extraction with a deepwell microplate.

7.4 DNA Extraction – tube format

- 7.4.1 Place lysis reagent (A+F) on a magnetic stir plate and stir on medium speed.

Note: A medium speed setting is used to ensure the lysis beads remain in suspension while dispensing. Using non-homogenous reagent could impact final results
--

- 7.4.2 Add 100 µl of the lysis reagent to each 1.5 ml screw cap conical tube.
- 7.4.3 Collect 100 µl of enriched sample and add to a tube containing lysis reagent. Mix the solution by pipetting up and down in the tube a minimum of 3 times.
- 7.4.4 Place tube in the heat block at 95-100°C for 15 min.
- 7.4.5 Vortex at high speed.
- 7.4.6 Centrifuge at 10,000-12,000 × g for 2 min.
- 7.4.7 Bring samples forward to 7.6.3.

7.5 DNA Extraction – deepwell format

- 7.5.1 Add 100 µl of the lysis reagent (A+F) to each well of a deepwell microplate.
- 7.5.2 After sample enrichment, collect a 100 µl aliquot of each sample and add into the

deepwell microplate.

- 7.5.3 For running iQ-Check™ Prep automation system, place unsealed deepwell plate, reagents and consumables into the instrument for automated sample processing.
- 7.5.4 For manual deepwell extraction, seal plate with pre-pierced sealing film.
- 7.5.5 Place in the Thermomixer® instrument at 95-100°C at 1,300 rpm for 15-20 min.
- 7.5.6 Remove plate from Thermomixer® and allow to cool in the refrigerator for a minimum of 10 minutes.
- 7.5.7 Bring samples forward to 7.6.3

7.6 PCR

- 7.6.1 Prepare PCR mix combining the amplification solution (Reagent C) and the fluorescent probes (Reagent B) according to the chart in the package insert.
- 7.6.2 Add 45 µl of PCR mix to each well of a 48 or 96 well PCR plate, pipetting carefully to avoid bubbles.
- 7.6.3 Add 5 µl of each sample, or control, to the appropriate well. If the DNA extract was obtained using the tube format (7.4.6), carefully pipette the supernatant, being careful not to disturb the pellet.
- 7.6.4 In addition to the laboratory's standard broth/enrichment controls, one positive control (Reagent E) and one negative control (Reagent D) should be run with each assay. Controls are provided as ready-to-use solutions in the kit.

Note: Steps 7.6.1 – 7.6.4 are completed automatically by the iQ-Check™ Prep automation system.

- 7.6.5 Seal the plate with the optical caps or use the PX1 plate sealer (for hard shell plates only).
- 7.6.6 Place the plate into the reaction module of the thermal cycler and close the lid.
- 7.6.7 Verify the plate setup and protocol selection is correct.
- 7.6.8 Click on the “run” button in the software program.
- 7.6.9 A window will appear allowing you to save your data in a specific location. Once “save” is selected, the PCR protocol is started.
- 7.6.10 A report will be generated at the end of the run assigning each sample a positive or negative result.

7.7 Interpretation of Results

- 7.7.1 CFX Manager IDE software will make the determination if a sample is positive, negative, inhibited or invalid.
- 7.7.2 A positive *Listeria monocytogenes* sample must have a Cq value ≥ 10 for the FAM fluorophore.

- 7.7.3 If there is no Cq value (Cq=N/A) for FAM, the internal control for that sample must then be analyzed. This sample is considered as a negative *Listeria monocytogenes* sample if the internal control has a Cq \geq 28.
- 7.7.4 If the internal control does not have a Cq value (Cq= N/A), then it is not possible to interpret the result. Such a result indicates an inhibition of the PCR reaction.
- 7.7.5 If a sample is scored as invalid, PCR should be repeated. If the sample is invalid a second time, contact Bio-Rad technical support and proceed with the confirmation steps for this sample.

7.8 Confirmation of Presumptive Positive Results

Using the 24 hour LSB enrichment broth aliquot, proceed with the plating and confirmation steps as described in the “isolation procedure” section of the cultural method MFHPB-30. Use at least one chromogenic agar listed in the cultural method in addition to Oxford agar. Additionally, transfer 0.1 mL of the presumptive positive broth to modified Fraser Broth and proceed with incubation and plating from Fraser Broth as indicated in MFHPB-30. Place the remainder of the 24 h LSB enriched sample back in the incubator and continue incubation for up to 48 h as required, for use as required in the confirmation of presumptive positive results, as described in MFHPB-30.

Once the presumptive positive has been confirmed, other confirmation streams can be discontinued. However, all plating and confirmation steps must be performed as described above, before reporting a presumptive positive result as a negative.

8. References

- 8.1 Bio-Rad Laboratories – iQ-Check™ *Listeria monocytogenes* II kit user guide. Available at <http://www.bio-rad.com/en-us/product/iq-check-listeria-monocytogenes-ii-kit>
- 8.2 Pagotto, F., Hébert, K., J. Farber. 2011. Isolation of *Listeria monocytogenes* and other *Listeria* spp. from foods and environmental samples (MFHPB-30). Volume 2. The *Compendium of Analytical Methods*.
- 8.3 Warburton, D., Boville, A., Pagotto, F., Daley, E. and Chow, C. 2012. The isolation of *Listeria monocytogenes* and other *Listeria* spp. from foods and environmental samples using Palcam broth (MFHPB-07). Volume 2. The *Compendium of Analytical Methods*

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