Health Products and Food Branch

Ottawa

**Detection of *Salmonella* spp. from all foods and selected environmental surfaces using
iQ-Check™ *Salmonella* Real-Time PCR Test Kit**

**Microbiological Methods Committee
Microbiology Evaluation Division
Bureau of Microbial Hazards, Food Directorate,
Health Products and Food Branch, Health Canada
Postal Locator: 2204E
Ottawa, Ontario K1A 0K9**

E-mail micro_methods_committee@hc-sc.gc.ca

1. Application

This method is applicable to the rapid detection of *Salmonella* spp. to determine compliance with the requirements of Sections 4 and 7 of the *Food and Drugs Act*. Positive results must be confirmed with a cultural method. This method has been validated for use in all foods, and the following environmental samples: stainless steel, plastic, ceramic and concrete.

2. Description

The iQ-Check™ system is a real-time PCR based assay that involves simultaneous amplification and detection of target organism after a single 21 ± 1 hr enrichment in Buffered Peptone Water. A regrowth step is not required. Results are automatically determined by the iQ-Check™ software. The assay should only be used by trained laboratory personnel who follow standard good laboratory practices.

3. Principle

The iQ-Check™ *Salmonella* kit is based on gene amplification and detection by real-time PCR. Ready-to-use PCR reagents contain DNA polymerase, nucleotides, DNA primers and a fluorescent DNA hybridization probe specific to a highly conserved region on the genome of all species of *Salmonella*. All reagents necessary for the test are contained in the kit. PCR is a technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling allow DNA denaturation, by heat, followed by primers binding to the target region. The DNA polymerase then uses these primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons. In real-time PCR, specific oligonucleotide probes are used to detect the DNA during the amplification, by hybridizing to the amplicons. These probes are linked to a fluorophore which fluoresces only when hybridized to the target sequence. In the iQ-Check™ *Salmonella* kit, carboxyfluorescein (FAM) is the fluorophore linked to the probe hybridizing to the *Salmonella* specific DNA sequence. In the absence of target DNA, no fluorescence will be detected, and the

sample is determined to be negative by the iQ-Check™ software. As the amount of amplicons increases with each round of amplification, fluorescence intensity also increases. During each PCR cycle, at the annealing step, the real-time PCR system measures this fluorescence and the associated software plots the fluorescence intensity versus number of cycles. This method allows a simple determination of the presence of *Salmonella* in a sample. To monitor for a successful DNA amplification in each reaction well, a synthetic DNA "internal control" is included in the reaction mix. This control is amplified with a specific probe at the same time as the *Salmonella* target DNA sequence, and detected by a second fluorophore.

4. Definition of Terms

See Appendix A of Volume 3.

5. Collection of Samples

See Appendix B of Volume 3

6. Materials and Special Equipment

Note: The Laboratory Supervisor must ensure that the analysis described in this method is carried out in accordance with the International Standard referred to as "ISO/IEC 17025:2005 (or latest version): General Requirements for the Competence of Testing and Calibration Laboratories".

Note: If the laboratory uses any variations of the media listed here (either product that is commercially available or made from scratch), it is their responsibility to ensure equivalency. Please forward equivalency data to the Editor of *Compendium of Analytical Methods* for consideration of modification of this method.

1) iQ-Check™ *Salmonella* kit (Bio-Rad catalog # 357-8123)

Reagent A - Lysis reagent, 1 bottle (20 mL), contains magnetic stirbar

Reagent B - Fluorescent probes, 1 tube (0.55 mL)

Reagent C - Amplification mix, 2 tubes (2.2 mL each)

Reagent D - PCR negative control, 1 tube (0.5 mL)

Reagent E - PCR positive control, 1 tube (0.25 mL)

2) Supplies and reagents

2.1 Buffered Peptone Water (BPW)

2.2 Conical screw cap sterile tubes, 1.5 mL (for tube extraction only)

2.3 Dey-Engley (D/E) neutralizing broth or Lethen broth for sponges and swabs

2.4 12.5 cm x 8.5 cm x 4 cm Deepwell microplates (for deepwell extraction only)

2.5 Low-profile unskirted PCR plates, 48 or 96 well, white (Bio-Rad catalog# MLL-4851XTU, MLL-9651XTU)

2.6 Optical Flat 8-cap strips, ultraclear (Bio-Rad catalog # TCS-0803XTU)

2.7 Pre-pierced sealing film (Bio-Rad catalog # 360-0040) or equivalent (for deepwell extraction only)

2.8 Powder-free gloves

2.9 Sterile filter tips - adaptable to 20 µL, 200 µL and 1000 µL micropipettes

2.10 Stomacher bag with incorporated filter

2.11 Micropipettes—20 µL, 200 µL and 1000 µL

3) Equipment

- 3.1 Bench top centrifuge - maximum 10,000-12,000 × g, for 1.5 mL tubes (for tube extraction only)
- 3.2 CFX Manager Software Industrial Diagnostic Edition (Bio-Rad catalog # 360-0042)
- 3.3 Dry heat block, for 1.5 mL tubes - capable of maintaining 95 - 100°C (for tube extraction only)
- 3.4 Incubator - capable of maintaining temperatures of 37 ± 1°C for sample incubation
- 3.5 Magnetic stir plate
- 3.6 Real-time PCR system – CFX96™ real-time PCR detection system (Bio-Rad catalog # 360-0037) or MiniOpticon™ real-time PCR detection system (Bio-Rad catalog # 359-1592)
- 3.7 Stomacher, blender or equivalent
- 3.8 Thermomixer (Bio-Rad catalog # 360-0038 [base] and 360-0039 [deepwell block]) or equivalent (for deepwell extraction only)
- 3.9 Vortex apparatus

Note:	It is the responsibility of each laboratory to ensure that the temperatures of the incubators or water baths are maintained at the recommended temperatures. Where 35°C is recommended in text of the method the incubator may be at 35 ± 1.0°C. Similarly, lower temperatures of 30 or 25 may be ± 1.0°C. However, where higher temperatures are recommended, such as 43 or 45.5°C, it is imperative that the incubators or water baths be maintained within 0.5°C due to potential lethality of the higher temperatures on the microorganism(s) being isolated.
--------------	---

7. Procedure

The test shall be carried out in accordance with the following instructions:

7.1 Handling of Sample Units

Note:	It is important that the D/E neutralizing or Lethen broth is used to moisten sponges and swabs. The use of other neutralizing buffers can cause PCR inhibition.
--------------	---

- 7.1.1 Add D/E neutralizing broth to sterile sponges and swabs to moisten.
- 7.1.2 The protocol of MFLP-41 should be followed for obtaining and transporting swabs.
- 7.1.3 In the laboratory prior to analysis, except for shelf-stable foods, keep food sample units refrigerated or frozen, depending on the nature of the product. Thaw frozen food samples in a refrigerator, or under time and temperature conditions which prevent microbial growth or death.
- 7.1.4 Analyze the sample units as soon as possible after receipt at the laboratory.

7.2 Preparation for Analysis

- 7.2.1 Turn on heating block or thermomixer in advance to reach 95 - 100°C.
- 7.2.2 Turn on the thermal cycler and computer.

7.3 Preparation of Samples

- 7.3.1 Temper BPW to approximately 37°C before use.
- 7.3.2 Weigh 25g or 25ml of each food sample and combine with 225 ml BPW. Blend, stomach or vortex as required for thorough mixing.
- 7.3.3 Add the environmental sponge or large swabs to 100 mL of BPW or composite up to 10 sponges with 100 mL of BPW per sponge. Place smaller environmental swabs (e.g., cotton tip) in 10 mL of BPW or composite up to 10 swabs with 10 mL BPW per swab. Sponges should be stomached for 2 min on high speed; swabs should be vortexed for 2 min on high speed.
- 7.3.4 Incubate all samples 37°C for 21 ± 1 hr.

7.4 DNA Extraction – tube format

- 7.4.1 Place lysis reagent on a magnetic stir plate and stir on medium speed.
- 7.4.2 Add 100 µl of the lysis reagent to each 1.5 ml screw cap conical tube.
- 7.4.3 Collect 100 µl of enriched sample and add to a tube containing lysis reagent. Mix the solution by pipetting up and down in the tube a minimum of 3 times.
- 7.4.4 Place tube in the heat block at 100°C for 15 min.
- 7.4.5 Vortex at high speed.
- 7.4.6 Centrifuge at 10,000-12,000 × g for 2 min.

7.5 DNA Extraction – deepwell format

- 7.5.1 Add 100 µl of the lysis reagent to each well of a deepwell microplate.
- 7.5.2 After sample enrichment, collect a 100 µl aliquot of each sample and add into the deepwell microplate.
- 7.5.3 Seal with pre-pierced sealing film.
- 7.5.4 Place in the thermomixer instrument at 99°C at 1,300 rpm for 10-15 min.
- 7.5.5 Remove plate from thermomixer and allow to cool in the refrigerator for 10 minutes.

7.6 PCR

- 7.6.1 Prepare PCR mix combining the amplification solution (Reagent C) and the fluorescent probes (Reagent B) according to the chart in the package insert.
- 7.6.2 Add 45 µl of PCR mix to each well of a 48 or 96 well PCR plate, pipetting carefully to avoid bubbles.
- 7.6.3 Add 5 µl of each sample, or control, to the appropriate well.
- 7.6.4 One positive control (Reagent E) and one negative control (Reagent D) should be run with each assay. Controls are provided as ready-to-use solutions in the kit.
- 7.6.5 Seal the plate with the optical caps.
- 7.6.6 Place the plate into the reaction module of the thermal cycler and close the lid.
- 7.6.7 Verify the plate setup and protocol selection is correct.
- 7.6.8 Click on the “run” button in the software program.
- 7.6.9 A window will appear allowing you to save your data in a specific location. Once “save” is selected, the PCR protocol is started.
- 7.6.10 A report will be generated at the end of the run assigning each sample a positive or negative result.

7.7 Confirmation of Presumptive Positive Results

Using the enrichment broth, immediately proceed with MFHPB-20, beginning with inoculation of the selective enrichment (Tetrathionate Brilliant Green broth and Rappaport-Vassiliadis Soya Peptone broth).

8. References

- 8.1 Health Canada, 2009. MFHPB-20. Methods for Isolation and Identification of *Salmonella* from Foods and Environmental Samples. Compendium of Analytical Methods, Volume 2. Available at <http://www.hc-sc.gc.ca/fn-an/res-rech/analy-meth/microbio/volume2/mfhp20-01-eng.php>

9. Preparation of Media

9.1 Letheen Broth

Beef Extract	5.0 g
Proteose Peptone no. 3	10.0 g
Polysorbate 80 (Tween™ 80)	5.0 g
Lecithin	0.7 g
Sodium chloride	5.0 g
Distilled water	1.0 L

Dissolve ingredients with gentle heating, if necessary, and adjust pH to 7.0 ± 0.2 . Dispense and autoclave for 15 min at 121°C.

9.2 Dey-Engley Neutralizing Broth

Tryptone	5.0 g
Yeast Extract	2.5 g
Dextrose	10.0 g
Sodium thioglycollate	1.0 g
Sodium thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$)	6.0 g
Sodium bisulfate (NaHSO_3)	2.5 g
Polysorbate 80 (Tween™ 80)	5.0 g
Lecithin (soy bean)	7.0 g
Bromcresol Purple	0.02 g
Distilled water	1.0 L

Dissolve ingredients with gentle heating, if necessary, and adjust pH to 7.6 ± 0.2 . Dispense and autoclave for 15 min at 121°C.

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